



19706 77<sup>th</sup> Avenue East, Bradenton, FL 34202, U.S.A.  
(941) 322-9739, DEW@skeeta.com, [www.skeeta.biz](http://www.skeeta.biz), CAGE: 3E1W4

## MEDICAL SHELTER (PN131-300)

### ASSEMBLY SEQUENCE:



[Medical Shelter packaged in pouch]



1. Remove the poles and reverse the pouch inside out to remove the Medical Shelter.



2. Construct frame by connecting (2) 60" & (2) 80" horizontal poles together using corner angles and position rectangular frame on floor, then insert (4) 80" vertical poles in lower corners, finally attached (2) 60" & (2) 80" horizontal poles together at top corner connectors.



3. Drape Medical Shelter over the top of the frame, while positioning the upper corners of Medical Shelter over the upper corners of the frame.



4. Temporarily disconnect 60" pole on the lower frame and slide pole through the lower end sleeve of Medical Shelter, and then reconnect pole at the lower corners.



5. Repeat step 4 on other end of the lower rectangular frame, to complete assembly.



6. Medical Shelter has zippered door flaps that wrap around the vertical poles to be binder clipped open.

### REPACKAGING SEQUENCE:

1. Unzip, to loosen Medical Shelter material and temporarily disconnect 60" horizontal pole sections from corner angles of lower frame to slide poles out of the sleeves.
2. Slide Medical Shelter up and over, away from the top of the frame to remove.
3. Disconnect 80" horizontal pole sections from corner angles of upper frame, while leaving angle connectors attached to each end of 60" pole sections. Then repeat disassembly until all the poles are loose.



4. Reverse the Medical Shelter inside out, zip closed and lay flat on ground with 80" pole sections on top at other end than pouch.



5. Roll-up Medical Shelter around 80" poles and position above pouch, then reverse pouch around, while stuffing Medical Shelter inside. Insert 60" poles with connectors on top, then zip pouch closed.